



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Platelet Receptor Gi24/VISTA/B7-H5/PD-1H Protein(C-Fc)

产品货号: TP04674

第三版

别名:Platelet receptor Gi24;Stress-induced secreted protein-1;Sisp-1;C10orf54;SISP1

描述:platelet receptor Gi24 is a single-pass type I membrane protein, and located at the cell surface. The protein can be cleaved by MMP14, and stimulate MMP14-mediated MMP2 activation. It is participated in the BMP signaling pathway. It also regulates the CD4-positive, alpha-beta T cell proliferation, and T cell cytokine production negatively. However, the protein can regulate stem cell differentiation positively.

配方:PBS

分子量:45.3 kDa

序列:Phe33-Ala194

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.