

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ubiquitin-Conjugating Enzyme E2 R2/UBE2R2/UBC3B Protein(N-6His)

产品货号: TP04675

第三版

别名: Ubiquitin-Conjugating Enzyme E2 R2; Ubiquitin Carrier Protein R2; Ubiquitin-Conjugating Enzyme E2-CDC34B; Ubiquitin-Protein Ligase R2; UBE2R2; CDC34B; UBC3B

描述: Ubiquitin-Conjugating Enzyme E2 R2 (UBE2R2) is a modification enzyme that belongs to the ubiquitin-conjugating enzyme family. UBE2R2 is involved in cell growth and transformation. It accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro, UBE2R2 catalyzes monoubiquitination and 'Lys-48'-linked polyubiquitination. It may be involved in degradation of katenin.

配方: PBS

分子量: 29.3 kDa

序列: Met 1-Ser238

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.