

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human L-selectin /sell Protein

产品货号: TP04681

第三版

别名:L-selectin (SELL); CD62 antigen-like family member L; Leukocyte adhesion molecule 1; Leukocyte surface antigen Leu-8; Leukocyte-endothelial cell adhesion molecule 1; Lymph node homing receptor; TQ1; gp90-MEL; LNHR and LYAM1

描述:SELL contains one C-type lectin domain, one EGF-like domain and two Sushi (CCP/SCR) domains. SELL is expressed in B-cell lines and T-lymphocytes. SELL functions as a cell surface adhesion protein and mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes. In addition, SELL also promotes initial tethering and rolling of leukocytes in endothelia.

配方:PBS

分子量:34.0 kDa

序列:Trp39-Asn332

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.