

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Transferrin Protein

产品货号: TP04687

第三版

别名: Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin; TF

描述: Serotransferrin belongs to transferrin family, and contains 2 transferrin-like domains. The protein is a secreted protein, and expressed by the liver and secreted in plasma. Transferrins are iron binding transport proteins which can bind two Fe³⁺ ions in association with the binding of an anion. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization. Serum transferrin may also have a further role in stimulating cell proliferation.

配方: PBS

分子量: 76.2 kDa

序列: Val20-Pro698

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.