

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interleukin-10 receptor subunit beta Protein

产品货号: TP04688

第三版

别名: Interleukin-10 receptor subunit beta(IL10RB); Cytokine receptor class-II member 4; Cytokine receptor family 2 member 4; Interleukin-10 receptor subunit 2

描述: IL10RB is a single-pass type I membrane protein and contains two fibronectin type-III domains. It is an accessory chain which is essential for the active interleukin 10 receptor complex. Coexpression of IL10RB and IL10RA proteins has been shown to be required for IL10-induced signal transduction. Defects in IL10RB are the cause of inflammatory bowel disease type 25 (IBD25) which is a chronic, relapsing inflammation of the gastrointestinal tract with a complex etiology.

配方: PBS

分子量: 24.5 kDa

序列: Met20-Ser220

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.