

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interleukin-13 Receptor Subunit Alpha-1 Protein

产品货号: TP04689

第三版

别名: Interleukin-13 receptor subunit alpha-1; IL-13 receptor subunit alpha-1; IL-13R subunit alpha-1; IL-13R-alpha-1; IL-13RA1; Cancer/testis antigen 19; CT19; CD213a1; IL13RA1; IL13R; IL13RA

描述: Interleukin-13 receptor subunit alpha-1 (IL13RA1) is a subunit of the interleukin 13 receptor. This subunit forms a receptor complex with IL4 receptor alpha, a subunit shared by IL13 and IL4 receptors. The human IL13-R 1 was originally cloned based on sequence homology to the mouse IL13-R 1, it share 76% aa sequence identity. Human The IL13-R 1 cDNA encodes a 427 amino acid (aa) residue precursor protein with a putative 21 aa residue signal peptide, a 324 aa residue extracellular domain, a 23 aa residue transmembrane region and a 59 aa residue cytoplasmic tail. The extracellular domain of IL13-R 1 is also closely related to that of IL13-R 2. It binds with low affinity to interleukin-13(IL13). IL13RA1 serves as a primary IL13- binding subunit of the IL13 receptor, and may also be a component of IL4 receptors. This protein has been shown to bind tyrosine kinase TYK2, and thus may mediate the signaling processes that lead to the activation of JAK1, STAT3 and STAT6 induced by IL13 and IL4.

配方: PBS

分子量: 37.7 kDa

序列: Gly22- Thr343

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.