

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ferritin Heavy Chain/FTH Protein(N-6His)

产品货号: TP04699

第三版

别名: Ferritin heavy chain; FTH1; FTH; FTHL6; Ferritin H subunit; Cell proliferation-inducing gene 15 protein

描述: Ferritin heavy polypeptide 1 (FTH1), is a ubiquitous intracellular protein which stores iron in a soluble, non-toxic, readily available form. FTH1 has ferroxidase activity and is important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Ferritin is composed of 24 subunits of the light and heavy ferritin chains. It plays a role in delivery of iron to cells and mediates iron uptake in capsule cells of the developing kidney. Variation of ferritin subunit composition may affect iron absorption and release in different tissues. Deficiency of ferritin proteins may cause several neurodegenerative diseases. Almost all living organisms can produce this protein, including algae, bacteria, higher plants, and animals.

配方: PBS

分子量: 23.4 kDa

序列: Met 1-Ser183

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.