



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Low-Density Lipoprotein Receptor Related Protein 12/LRP12 Protein(C-6His)

产品货号: TP04700

第三版

别名:Low-Density Lipoprotein Receptor-Related Protein 12; LRP-12; Suppressor of Tumorigenicity 7 Protein; LRP12; ST7

描述:Low-Density Lipoprotein Receptor-Related Protein 12 (LRP12) belongs to the LDLR family. LRP12 is a type I transmembrane protein and widely expressed in heart, skeletal muscle, brain, lung, placenta and pancreas. LRP12 contains 2 CUB domain and 5 LDL-receptor class A domain. LRP12 has been shown to interact with GNB2L1, ZFYVE9 and ITGB1BP3. LRP12 is a receptor probably, which may be involved in the internalization of lipophilic molecules and/or signal transduction. In addition, LRP12 may act as a tumor suppressor.

配方:PBS

分子量:52.5 kDa

序列:Asn28-Ile488

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.