

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Mammalian Ependymin-Related Protein 1/EPDR1/UCC1 Protein(C-6His)

产品货号: TP04713

第三版

别名:Mammalian ependymin-related protein 1;EPDR1;Upregulated in colorectal cancer gene 1 protein;MERP1

描述:EPDR1 is a member of the ependymin family. EPDR1 is a type II transmembrane protein that is similar to two families of cell adhesion molecules, the protocadherins and ependymins. It may play a role in calcium-dependent cell adhesion. EPDR1 is glycosylated, and the orthologous mouse protein is localized to the lysosome. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8.

配方:PBS

分子量:22.6 kDa

序列:Ala38-Ser223

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.