



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Microfibril-Associated Glycoprotein 4/MFAP4 Protein(C-6His)

产品货号: TP04718

第三版

别名: Microfibril-associated glycoprotein 4

描述: Microfibril-associated glycoprotein 4(MFAP4) is a secreted protein and contains 1 fibrinogen C-terminal domain, similarity to a bovine microfibril-associated protein. The protein has binding specificities for both collagen and carbohydrate. It is thought to be an extracellular matrix protein which is involved in calcium-dependent cell adhesion or intercellular interactions. The gene is located within the Smith-Magenis syndrome region.

配方: PBS

分子量: 27.5 kDa

序列: Val22-Ala255

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.