



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Chymotrypsinogen B/CTRB1 Protein(C-6His)

产品货号: TP04763

第三版

别名:Chymotrypsinogen B; CTRB1; CTRB

描述:Chymotrypsinogen B (CTRB1) is a 263 amino acid protein with signal peptide (1-18) and Chymotrypsinogen B (19-263). Chymotrypsinogen B have three chains:Chymotrypsin B chain A(19-31), Chymotrypsin B chain B(34-164), Chymotrypsin B chain C (167-263). Chymotrypsinogen B is a Serine Protease Hydrolase secrets into gastrointestinal tract as the inactive precursor Chymotrypsinogen.

配方:PBS

分子量:26.9 kDa

序列:Cys19-Asn263

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.