

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Xaa-Pro Aminopeptidase 2 Protein

产品货号: TP04771

第三版

别名:Xaa-Pro Aminopeptidase 2; Aminoacylproline Aminopeptidase; Membrane-Bound Aminopeptidase P; Membrane-Bound APP; Membrane-Bound AmP; mAmP; X-Pro Aminopeptidase 2; XPNPEP2

描述:Xaa-Pro aminopeptidase 2 (XPNPEP2) belongs to the peptidase M24B family of metalloproteases. Human XPNPEP2 is widely expressed in various tissues, such as kidney, lung, heart, placenta, liver, small intestine and colon. However, it doesn't express in brain, skeletal muscle, pancreas, spleen, thymus, prostate, testis and ovary. XPNPEP2 is a Homotrimer which binds 2 manganese ions per subunit. The metalloprotease XPNPEP2 may play a role in the inflammatory process and other reactions produced in response to injury or infection and the metabolism of the vasodilator bradykinin.

配方:PBS

分子量:71.6 kDa

序列:His22-Ala650

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.