



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein FAM3C Protein(C-6His)

产品货号: TP04779

第三版

别名:Protein FAM3C; Interleukin-Like EMT Inducer; FAM3C; ILEI

描述:FAM3C, also called interleukin-like EMT inducer, usually exist in most secretory epithelia. It belongs to the FAM3 family according to their sequence similarities. The up-regulation and/or mislocalization in breast cancer and liver carcinoma cells of FAM3C is strongly correlated with metastasis formation and survival. FAM3C can be involved in retinal laminar formation and promote epithelial to mesenchymal transition.

配方:PBS

分子量:23.2 kDa

序列:Gln25-Asp227

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.