
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse / Recombinant Human / Rat / Cynomolgus / Canine BDNF Protein

产品货号: TP04789

第三版

别名:BDNF

描述:BDNF is a member of the nerve growth factor family. It is highly expressed in hippocampus, amygdala, cerebral cortex and cerebellum. It also can be detected in heart, lung, skeletal muscle, testis, prostate and placenta. BDNF is induced by cortical neurons, and is necessary for survival of striatal neurons in the brain. During development, BDNF promotes the survival and differentiation of selected neuronal populations of the peripheral and central nervous systems. It participates in axonal growth, pathfinding and in the modulation of dendritic growth and morphology. It functions as the major regulator of synaptic transmission and plasticity at adult synapses in many regions of the CNS. The versatility of BDNF is emphasized by its contribution to a range of adaptive neuronal responses including long-term potentiation (LTP), long-term depression (LTD), certain forms of short-term synaptic plasticity, as well as homeostatic regulation of intrinsic neuronal excitability.

配方:PBS

分子量:13.5 kDa

序列:Met 1-Arg249

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.