
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Uteroglobin / SCGB1A1 Protein (His tag)

产品货号: TP04801

第三版

别名:CC10;CC16;CCSP;PCB-BP;UG;UGB;Utg

描述:Uteroglobin (UG), also known as Secretoglobin 1A member 1 (SCGB1A1), Blastokinin, Clara cell secretor protein (CCSP) or Clara cell-specific 10-kDa protein (CC10), is the founding member of the secretoglobin family of small, secreted, disulfide-bridged dimeric proteins found only in mammals. This protein is mainly expressed in lung, with anti-inflammatory/immunomodulatory properties. Previous in vitro studies demonstrated that CCAAT/enhancer-binding proteins (C/EBPs) are the major transcription factors for the regulation of SCGB1A1 gene expression, whereas FOXA1 had a minimum effect on the transcription. Uteroglobin is a multifunctional protein with antiinflammatory/immunomodulatory properties. Uteroglobin inhibits soluble phospholipase A(2) activity and binds and perhaps sequesters hydrophobic ligands such as progesterone, retinols, polychlorinated biphenyls, phospholipids, and prostaglandins. In addition to its antiinflammatory activities, Uteroglobin manifests antichemotactic, antiallergic, antitumorigenic, and embryonic growth-stimulatory activities. The tissue-specific expression of the Uteroglobin gene is regulated by several steroid hormones, although a nonsteroid hormone, prolactin, further augments its expression in the uterus. Based on its anti-inflammatory and antiallergic properties, Uteroglobin is a potential drug target. The mechanism of Uteroglobin action is likely to be even more complex as it also functions via a putative receptor-mediated pathway.

配方:PBS

分子量:9.8 kDa

序列:Met 1-Phe 96

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.