

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human NKG2D Ligand 1/NKG2DL/ULBP1 Protein(C-Fc)

产品货号: TP04819

第三版

别名:NKG2D ligand 1;NKG2DL1;ALCAN-beta;Retinoic acid early transcript 1I;UL16-binding protein 1;ULBP1

**描述:**NKG2D ligand 1, also called ULBP1, is a member of UL16-binding protein (ULBP) family which has also been termed the retinoic acid early transcript 1 (RAET1) family. Unlike the classical MHC class I molecules and the MIC molecules possess 1, 2 and 3 domains, ULBP/RAET1 family members lack 3 domain. ULBP1 is recognized by the activating receptor NKG2D on the surface of cytotoxic natural killer (NK) and T cells, and then promotes the lysis of cells expressing ULBP1 which is important for the immune surveillance. ULBP1 and several other family members, ULBP2 and ULBP5, own the ability to bind the human cytomegalovirus (CMV) UL16 glycoprotein. The human CMV glycoprotein UL16 binds to intracellular ULBP1 and so inhibits its expression at the cell surface, which reduces the susceptibility of the virus-infected cell to cytotoxic destruction by NK cells. The expression of ULBP1 has been found on some tumor cells and is implicated in tumor surveillance.

**配方:**PBS

**分子量:**49.4 kDa

**序列:**Gly26-Pro215

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.