

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Follistatin-Like 1/FSTL1 Protein(C-Fc)

产品货号: TP04823

第三版

别名:Follistatin-related protein 1; Follistatin-like protein 1; TGF-beta-inducible protein TSC-36; Fstl1

描述:Follistatin-like 1 (FSTL1) is a secreted glycoprotein that has been grouped into the follistatin family of proteins. FSTL1 is composed of a follistatin domain and two non-functional calcium-binding motifs. It was originally cloned as a TGF- β 1 inducible factor but subsequently shown to regulate diverse developmental pathways and tissue homeostasis. Ablation of the FSTL1 gene in the mouse results in several structural developmental defects and neonatal lethality due to respiratory failure. FSTL1 suppresses BMP signaling, but the precise mechanism of its action has not been elucidated. FSTL1 is expressed in the human placenta, mainly in extravillous trophoblasts.

配方:PBS

分子量:59.6 kDa

序列:Glu19-Ile306

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.