

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Frizzled 1/FZD1 Protein(C-Fc)

产品货号: TP04824

第三版

别名:Frizzled-1; Fz-1; mFz1; Fzd1; Frizzled homolog 1

描述:Frizzled homolog 1 (FZD1), also known as Frizzled-1, a member of the G-protein-coupled receptor superfamily, mediates Wnt signaling by serving as a co-receptor with LRP-5 for Wnt ligands. The FZD1 protein contains a signal peptide, a cysteine-rich domain in the N-terminal extracellular region, 7 transmembrane domains, and a C-terminal PDZ domain-binding motif.

配方:PBS

分子量:46.8 kDa

序列:Val69-His248

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.