



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Tissue Factor Pathway Inhibitor 2/TFPI-2/PP5/REF-1 Protein(C-6His)

产品货号: TP04857

第三版

别名:Tissue Factor Pathway Inhibitor 2; TFPI-2; Placental Protein 5; PP5; TFPI2

描述:Human Tissue Factor Pathway Inhibitor 2 (TFPI2) has a N-terminal acidic region, three Kunita domains separated with by two linker regions, and a C-terminal basic region. TFPI2 has the function of regulating plasmin-mediated matrix remodeling, inhibits trypsin, plasmin, factor VIIa/tissue factor and weakly factor Xa.. TFPI2 has no effect on thrombin. TFPI2 may contribute to tumor progression in many cancers; in these cancers, the expression of TFPI2 can be down-regulated.

配方:PBS

分子量:22.9 kDa

序列:Asp23-Lys213

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.