

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cystatin A Protein(N-6His)

产品货号: TP04864

第三版

别名:Cystatin-A; Cystatin-AS; Stefin-A; CSTA; STF1; STFA

描述:Human Cystatin A (CSTA) is a member of family 1 of the cystatin superfamily, which is characterized by lacking of disulphide bonds and carbohydrates. Cystatin A is an intracellular inhibitor regulating the activities of cysteine proteases of the papain family such as Cathepsins B, H and L. Cystatin A is also implicated in a number of disease states. Due to altered proteolytic state in cancer progression, Cystatin A may play a role in the proteolytic pathways.

配方:PBS

分子量:11.8 kDa

序列:Ile2-Phe98

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.