



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cystatin SA/CST2 Protein(C-6His)

产品货号: TP04867

第三版

别名:Cystatin-SA; Cystatin-2; Cystatin-S5; CST2

描述:Cystatin SA is a member of family 2 of the Cystatin superfamily. Together with Cystatins S and SN, it is produced by the salivary gland and secreted largely in the submandibular/sublingual saliva. Cystatin SA inhibits Cathepsin L, but may not inhibit Cathepsins B and H. The functions of Cystatin SA are largely unknown. However, it has been suggested to play a role in the control of periodontal tissue destruction.

配方:PBS

分子量:15.4 kDa

序列:Trp21-Ala141

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.