

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Fc Receptor-Like Protein 1 Protein

产品货号: TP04874

第三版

别名:Fc Receptor-Like Protein 1; FcR-Like Protein 1; FcRL1; Fc Receptor Homolog 1; FcRH1; IFGP Family Protein 1; hIFGP1; Immune Receptor Translocation-Associated Protein 5; CD307a; FCRL1; FCRH1; IFGP1; IRTA5

描述:Fc Receptor-Like Protein 1 (FCRL1) is a single-pass type I membrane protein that may function as an activating coreceptor in B cells. FCRL1 is primarily expressed in secondary lymphoid tissues by mature subsets of B cells. It can be detected in the spleen, lymph node, heart, skeletal muscle, kidney, liver and placenta. It is specifically expressed by mature B lineage cells with higher expression at the protein level in naive B cells compared to memory B cells. FCRL1 contains three extracellular C2-like immunoglobulin domains, a transmembrane domain, and a cytoplasmic domain with two immunoreceptor-tyrosine activation motifs and may play a role in the regulation of cancer cell growth.

配方:PBS

分子量:32.2 kDa

序列:Ala17-Asn303

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.