

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Anterior Gradient Protein 3 Homolog/AG-3/BCMP11/AGR3 Protein(C-6His)

产品货号: TP04875

第三版

别名:Anterior Gradient Protein 3 Homolog; AG-3; AG3; hAG-3; Breast Cancer Membrane Protein 11; AGR3; BCMP11

描述:Anterior Gradient Protein 2(AG-2) and Anterior Gradient Protein 3 (AG-3) are human homologues of genes involved in differentiation, are associated with oestrogen receptor-positive breast tumours and interact with metastasis gene C4.4a and dystroglycan (hAG-3 protein). AG-3 could serve as a prognostic marker for survival in patients with low grade and high grade serous ovarian carcinomas.

配方:PBS

分子量:18.0 kDa

序列:Ile22-Leu166

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.