
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human T-Cell Surface Glycoprotein CD8 beta Chain Protein

产品货号: TP04883

第三版

别名:T-Cell Surface Glycoprotein CD8 Beta Chain; CD8b; CD8B; CD8B1

描述:T-Cell Surface Glycoprotein CD8 Chain (CD8 Antigen) is a cell surface glycoprotein found on most cytotoxic T lymphocytes that mediates efficient cell-cell interactions within the immune system. CD8 Antigen, acting as a coreceptor, and the T-cell receptor on the T lymphocyte recognize antigens displayed by an antigen presenting cell (APC) in the context of class I MHC molecules. The functional coreceptor is either a homodimer composed of two alpha chains, or a heterodimer composed of one alpha and one beta chain. Both alpha and beta chains share significant homology to immunoglobulin variable light chains. Multiple alternatively spliced transcript variants encoding distinct membrane associated or secreted isoforms have been described. A pseudogene, also located on chromosome 2, has been identified.

配方:PBS

分子量:17.8 kDa

序列:Leu22-Pro170

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.