

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Bcl-2-Like Protein 2/BCL2L2 Protein(C-6His)

产品货号: TP04885

第三版

别名:Bcl-2-Like Protein 2; Bcl2-L-2; Apoptosis Regulator Bcl-W; BCL2L2; BCLW; KIAA0271

描述:Bcl-2-like protein 2 (BCL2L2) belongs to the Bcl-2 family. BCL2L2 is highly expressed in the brain, spinal cord, testis, pancreas, heart, spleen, and mammary glands. BCL2L2 is a peripheral membrane protein containing three motifs, BH1, BH2 and BH4. The BH4 motif appears to be involved in the anti-apoptotic function. The BH1 and BH2 motifs form a hydrophobic groove which acts as a docking site for the BH3 domain of some pro-apoptotic proteins. BCL2L2 promotes cell survival and blocks dexamethasone-induced apoptosis. Furthermore, BCL2L2 mediates survival of postmitotic Sertoli cells by suppressing death-promoting activity of BAX.

配方:PBS

分子量:19.9 kDa

序列:Ala2-Thr172

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.