

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protocadherin-10/PCDH10 Protein(C-Fc)

产品货号: TP04887

第三版

别名:Protocadherin-10; PCDH10; KIAA1400

描述:Protocadherin-10 (PCDH10) is a single-pass type I membrane protein that contains six extracellular cadherin domains, one transmembrane domain and one cytoplasmic tail differing from those of the classical cadherins. As a potential calcium-dependent cell-adhesion neuronal receptor, it may play a role in the establishment and function of specific cell-cell connections in the brain. PCDH10 moderately expressed in all regions of the brain examined, as well as in testis and ovary, and low expression in all other tested tissues.

配方:PBS

分子量:102.4 kDa

序列:Gln19-Thr715

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.