

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Asialoglycoprotein Receptor 1 Protein

产品货号: TP04891

第三版

别名:Asialoglycoprotein Receptor 1; ASGP-R 1; ASGPR 1; C-Type Lectin Domain Family 4 Member H1; Hepatic Lectin H1; HL-1; ASGR1; CLEC4H1

描述:Asialoglycoprotein Receptor 1 (ASGPR1) is an endocytic recycling receptor, belongs to the long-form subfamily of the C-type/Ca²⁺-dependent lectin family. ASGPR consists of two noncovalently-linked subunits, ASGPR1 and ASGPR2. ASGPR1 mediates the endocytosis of plasma glycoproteins, recognizes terminal galactose and N-acetylgalactosamine units. When the ligand binds to ASGPR1, results in the complex is internalized and transported to a sorting organelle, then ASGPR1 and ligand can be disassociated, ASGPR1 returns to the cell membrane surface.

配方:PBS

分子量:27.4 kDa

序列:Gln62-Ile291

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.