

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CMRF35-Like Molecule 9 Protein

产品货号: TP04930

第三版

别名:CMRF35-Like Molecule 9; CLM-9; CD300 Antigen-Like Family Member G; Triggering Receptor Expressed on Myeloid Cells 4; TREM-4; CD300g; CD300LG; CLM9; TREM4

描述:CMRF35-Like Molecule 9 (CD300LG) is a single-pass type I membrane protein which belongs to the CD300 family. CD300LG has one Ig-like V-type domain which mediates binding to lymphocyte. CD300LG is highly expressed in heart, skeletal muscle and placenta. CD300LG acts as a receptor which may mediate L-selectin-dependent lymphocyte rollings. CD300LG also binds SELL in a calcium dependent manner and lymphocyte. CD300LG may play a important role in molecular traffic across the capillary endothelium.

配方:PBS

分子量:25.8 kDa

序列:Leu19-Arg247

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.