

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Insulin-Like Growth Factor II/IGF2 Protein

产品货号: TP04933

第三版

别名: Insulin-Like Growth Factor II; IGF-II; Somatomedin-A; IGF2; PP1446

描述: Insulin-Like Growth Factor II (IGF2) belongs to the insulin family of polypeptide growth factors that is involved in development and growth. Members of this family are structurally homologous to proinsulin, and share higher sequence identity. IGF2 is expressed only from the paternally inherited allele and believed to be secreted by the liver and to circulate in the blood. IGF2 possess growth-promoting activity and can stimulate the proliferation and survival of various cell types including muscle, bone, and cartilage tissue in vitro. IGF2 is influenced by placental lactogen and may play a role in fetal development.

配方: PBS

分子量: 8.9 kDa

序列: Ala25-Glu91

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.