

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Secreted and Transmembrane Protein 1 Protein

产品货号: TP04937

第三版

别名: Secreted and Transmembrane Protein 1; Protein K-12; SECTM1; K12

描述: Secreted and Transmembrane Protein 1 (SECTM1) is a transmembrane and secreted protein that belongs to the SECTM family. SECTM1 is expressed in a perinuclear Golgi-like pattern. It is detected at the highest levels in peripheral blood leukocytes and breast cancer cell lines. SECTM1 is considered to participate in thymocyte signaling and the hematopoietic/immune system processes. It is reported that SECTM1 is a broadly expressed, IFN- γ -inducible molecule, which functions as a potent costimulatory ligand for T cell activation and is synergistic with anti-CD28.

配方: PBS

分子量: 13.8 kDa

序列: Gln29-Gly145

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.