



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse TACI/TNFRSF13B/CD267 Protein(C-Fc)

产品货号: TP04945

第三版

别名:Tumor necrosis factor receptor superfamily member 13B;Transmembrane activator and CAML interactor;CD267; Tnfrsf13b

描述:Tumor necrosis factor receptor superfamily member 13B, also known as TNFRSF13B or more commonly as TACI (transmembrane activator and CAML interactor), is a transmembrane receptor protein found predominantly on the surface of B cells, which are an important part of the immune system.

配方:PBS

分子量:41.4 kDa

序列:Phe5-Thr129

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.