

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cytokine receptor common subunit gamma Protein

产品货号: TP04947

第三版

别名:Cytokine receptor common subunit gamma;IL2RG;Interleukin-2 receptor subunit gamma;IL-2 receptor subunit gamma;IL-2R subunit gamma;gammaC

描述:IL2RG contains one fibronectin type-III domain. IL2RG is an important signaling component of many interleukin receptors, including those of interleukin -2, -4, -7 and -21, and is thus referred to as the common gamma chain. IL2RG interacts with SHB upon interleukin stimulation and HTLV-1 accessory protein p12I. Defects in IL2RG are the cause of X-linked combined immunodeficiency (XCID) and severe combined immunodeficiency X-linked T-cell-negative /B-cell-positive / NK-cell-negative (XSCID).

配方:PBS

分子量:55.4 kDa

序列:Leu23-Ala262

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.