

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human FGF R5/FGFRL1 Protein(C-6His)

产品货号: TP04948

第三版

别名:Fibroblast Growth Factor Receptor-Like 1; FGF Receptor-Like Protein 1; FGF Homologous Factor Receptor; FGFR-Like Protein; Fibroblast Growth Factor Receptor 5; FGFR-5; FGFRL1; FGFR5; FHFR

描述:Fibroblast Growth Factor Receptor-Like 1 (FGFRL1) is a single-pass type I membrane protein that belongs to the FGF receptor family. The mature human FGFRL1 consists of a 354 amino acid extracellular domain (ECD) with 3 Ig-like C2-type domains, a 21 amino acid transmembrane segment, and a 134 amino acid cytoplasmic domain. FGFR1 expressed in various tissues, preferentially in cartilaginous tissues and pancreas. It highly expressed in the liver, kidney, heart, brain and skeletal muscle, weakly expressed in the lung, small intestine and spleen. FGFRL1 has a negative effect on cell proliferation.

配方:PBS

分子量:39.9 kDa

序列:Ala25-Pro378

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.