



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Myeloid Differentiation Protein 1/Lymphocyte Antigen 86 Protein

产品货号: TP04954

第三版

别名:Lymphocyte antigen 86;LY86;MD1;Protein MD-1

描述:LY86 is a protein that in humans is encoded by the LY86 gene, which May cooperate with CD180 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) and cytokine production. Important for efficient CD180 cell surface expression.

配方:PBS

分子量:42.8 kDa

序列:Gly21-Ser162

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.