

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human UBE2T Protein (His Tag)

产品货号: TP04981

第三版

别名:HSPC150;PIG50

描述:Ube2T is the E2 ubiquitin-conjugating enzyme of the Fanconi anemia DNA repair pathway and it is overexpressed in several cancers, representing an attractive target for the development of inhibitors. Notably, UBE2T locates at 1q32.1, and the gain of 1q is frequently observed in a variety of cancers. For instance, UBE2T serves an important role in the growth of bladder cancer cells, and may be considered as a potential biomarker and therapeutic target for bladder cancer.

配方:PBS

分子量:23.9 kDa

序列:Met 1-Val197

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.