

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human CALML3 Protein (His tag)

产品货号: TP04986

第三版

别名:CLP

**描述:** Calmodulin-like protein 3 (CALML3) is similar to that of authentic calmodulin and may actually compete with calmodulin by binding, with different affinity, to cellular substrates. Calmodulin-like protein 3 (CALML3) is a tumor-sensitive protein specifically expressed in normal epithelial cells but downregulated in tumorigenesis. Downregulation of the protein is an early event in breast cancer development. One of the most pressing questions raised by the discovery of CLP/CALML3 is that of its potential targets. Although it is 85% identical to human calmodulin, the distinct properties of CLP suggest that it has specific targets or targets that only partially overlap with those of calmodulin. Research has identified the unconventional myosin-10 (Myo10) as a specific target of CALML3. The discovery of Myo10 as a specific target of CALML3 is highly significant and suggests multiple lines of further research such as investigations of the Ca<sup>2+</sup> regulation of Myo10 and the role of the loss of CLP in epithelial differentiation, adhesion, and cancer. Cells expressing CALML3 displayed a striking increase in the number and length of myosin-10-positive filopodia and showed increased mobility in a wound healing assay.

**配方:**PBS

**分子量:**18.7 kDa

**序列:**Met 1-Lys 149

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.