

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Growth Hormone/GH Protein

产品货号: TP04990

第三版

别名:GH1;Somatotropin;Growth hormone;GH;GH-N;Growth hormone 1;Pituitary growth hormone

**描述:**Growth hormone (GH), also known as somatotropin, is a member of a family of growth factors. It plays an important role in growth control. Its major role in stimulating body growth is to stimulate the liver and other tissues to secrete IGF-1. GH includes prolactin, placental lactogens, proliferins, and somatolactin. It is synthesized primarily by somatotropes in the anterior pituitary and is stored in secretory granules. It stimulates both the differentiation and proliferation of myoblasts. It also stimulates amino acid uptake and protein synthesis in muscle and other tissues.

**配方:**PBS

**分子量:**22.1 kDa

**序列:**Phe27-Phe217

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.