

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Neuritin 1-Like Protein/NRN1L Protein(C-6His)

产品货号: TP04998

第三版

别名:Neuritin-like protein; NRN1L; UNQ2446/PRO5725

**描述:**Neuritin-like protein belongs to the neuritin family. Neuritin is a GPI-anchored protein that promotes neurite outgrowth and branching of neuritic processes in primary hippocampal and cortical cells. Neuritin expression also enhances the development of motor neuron axon arbors by promoting neuromuscular synaptogenesis and by stimulating the addition of new axon branches. Neuritin is induced by neuronal activity and by the neurotrophins, BDNF and NT3. NRN1L contains a consensus cleavage signal found in glycosylphosphatidylinositol (GPI)-anchored proteins.

**配方:**PBS

**分子量:**12.3 kDa

**序列:**Ala36-Ala139

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.