



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Syntaxin-7/STX7 Protein(N-6His)

产品货号: TP05009

第三版

别名:Syntaxin-7; STX7; syntaxin 7

描述:Syntaxin 7 (STX7) is a member of a family of membrane-trafficking proteins named as soluble N-ethylmaleimide-sensitive factor attachment protein receptors (SNAREs), and mediates the endocytic trafficking from early endosomes to late endosomes and lysosomes. Syntaxin-7 has been shown to be present in both late endosomes and lysosomes, and to be required for both homotypic late endosome fusion and heterotypic fusion with lysosomes. STX7 has been shown to interact with STX8, VPS18, Vesicle-associated membrane protein 8 and VPS11.

配方:PBS

分子量:28.0 kDa

序列:Ser2-Leu238

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.