

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Syntaxin-8/STX8 Protein

产品货号: TP05010

第三版

别名:Syntaxin-8; STX8

**描述:**Syntaxin-8 is a single-pass type IV membrane protein which belongs to the syntaxin family. It contains one t-SNARE coil homology domain. STX8 is highly expressed in heart, also found in brain, kidney, liver, lung, placenta, skeletal muscle, spleen and pancreas. STX8 is involved in protein trafficking from early to late endosomes via vesicle fusion and exocytosis. It as a vesicle trafficking protein functions in the early secretory pathway, possibly mediating retrograde transport form cis-golgi membrane to the ER.

**配方:**PBS

**分子量:**24.6 kDa

**序列:**Met 1-Gly215

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.