



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human NKG2D Ligand 2/NKG2DL2/N2DL2 Protein(C-Fc)

产品货号: TP05024

第三版

别名:NKG2D Ligand 2; N2DL-2; NKG2DL2; ALCAN-Alpha; Retinoic Acid Early Transcript 1H; UL16-Binding Protein 2; ULBP2; N2DL2; RAET1H

描述:NKG2D Ligand 2 (N2DL2) is a member of a family of cell-surface proteins. N2DL2 function as ligands for human cytomegalovirus glycoprotein UL16. N2DL2 is anchored to the membrane via a GPI-linkage. N2DL2 is bind to human NKG2D, an activating receptor expressed on NK cells, NKT cells, T cells. Engagement of NKG2D results in the activation of cytolytic activity and cytokine production by these effects cells. The ULBPs are expressed on some tumor cells and have been implicated in tumor surveillance.

配方:PBS

分子量:51.4 kDa

序列:Gly26-Ser217

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.