

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse CD5 antigen-like Protein

产品货号: TP05027

第三版

别名:CD5 antigen-like;Apoptosis inhibitor expressed by macrophages;Apoptosis inhibitory 6;CT-2;SP-alpha;Cd5l;Aim; Api6

描述:CD5L, also known as CD5 antigen-like, is secreted protein which mainly expressed by macrophages in lymphoid and inflamed tissues and regulates mechanisms in inflammatory responses, such as infection or atherosclerosis. It is able to inhibit lipid droplet size in adipocytes. CD5L acts as a key regulator of metabolic switch in T-helper Th17 cells. It is Participates in obesity-associated autoimmunity via its association with IgM, interfering with the binding of IgM to Fc α / μ receptor and enhancing the development of long-lived plasma cells that produce high-affinity IgG autoantibodies. It also acts as an inhibitor of apoptosis in macrophages: promotes macrophage survival from the apoptotic effects of oxidized lipids in case of atherosclerosis. It is involved in early response to microbial infection against various pathogens by acting as a pattern recognition receptor and by promoting autophagy.

配方:PBS

分子量:37.4 kDa

序列:Glu22-VaL352

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.