
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human NPDC1/CAB1 Protein(C-6His)

产品货号: TP05029

第三版

别名:Neural proliferation differentiation and control protein 1;NPDC-1;NPDC1;RP11-229P13.1; CAB; CAB-1; CAB1

描述:Neural proliferation differentiation and control protein 1(NPDC1) is a protein that in humans is encoded by the NPDC1 gene. It is a single-pass membrane protein and belongs to the NPDC1/cab-1 family. The protein strongly expressed in adult brain and especially in hippocampus, frontal lobe and temporal lobe. The protein suppresses oncogenic transformation in neural and non-neural cells and down-regulates neural cell proliferation and it might be involved in transcriptional regulation.

配方:PBS

分子量:16.5 kDa

序列:Gly35-Asp181

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.