

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Copine-1/CPNE1 Protein(N, C-6His)

产品货号: TP05037

第三版

别名:Copine-1; Copine I; CPN1; CPNE1

描述:Copine-1(CPNE1) encodes a calcium-dependent protein which belongs to the copine family. CPNE1 contains two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus. However, CPNE1 does not contain a predicted signal sequence or transmembrane domains. CPNE1 may regulate molecular events at the interface of the cell membrane and cytoplasm. CPNE1 has a broad tissue distribution and it may function in membrane trafficking.

配方:PBS

分子量:62.2 kDa

序列:Met 1-Ala537

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.