

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Transcobalamin II Receptor/TCbIR/8D6A/CD320 Protein(C-6His)

产品货号: TP05042

第三版

别名:CD320 antigen;8D6 antigen;FDC-signaling molecule 8D6;FDC-SM-8D6;Transcobalamin receptor;TCbIR;CD320

描述:CD320 antigen is also known as 8D6 antigen , FDC-signaling molecule 8D6 , Transcobalamin receptor and 8D6A. It is a single-pass type I membrane protein and containing two LDL-receptor class A domains. CD320 has been recently discovered and reported as a follicular dendritic cell (FDC) protein. CD320 can augments the proliferation of plasma cells precursors generated by IL-10. CD320 also founctions a receptor for the cellular uptake of transcobalamin bound cobalamin. Defects in CD320 are the cause of methylmalonic aciduria type TCbIR (MMATC) which is a metabolic disorder.

配方:PBS

分子量:21.1 kDa

序列:Ser36-Val231

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.