

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Parathyroid Hormone 1 Receptor/PTH1R Protein(Gly49, C-6His)

产品货号: TP05048

第三版

别名:Parathyroid hormone 1 receptor;PTH1R

描述:Parathyroid hormone 1 receptor(PTH1R) is a multi-pass membrane protein. The protein is expressed in high levels in bone and kidney and regulates calcium ionhomeostasis through activation of adenylate cyclase and phospholipase C. In bone, it is expressed on the surface of osteoblasts. When the receptor is activated through PTH binding, osteoblasts express RANKL (Receptor Activator of Nuclear Factor kB Ligand), which binds to RANK (Receptor Activator of Nuclear Factor kB) on osteoclasts. This turns on osteoclasts to ultimately increase the resorption rate.

配方:PBS

分子量:20.2 kDa

序列:Tyr23-Met 189

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.