

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TREML2/TLT-2 Protein(C-Fc)

产品货号: TP05062

第三版

别名:Trem-like transcript 2 protein;TLT2; Triggering receptor expressed on myeloid cells-like protein 2; TLT2;C6orf76;

描述:Trem-like transcript 2 protein (TLT2), also known as Triggering receptor expressed on myeloid cells-like protein 2, TLT2 and C6orf76, is single-pass type I membrane protein. TREML2 contains one Ig-like V-type domain, which can be induced in CD4 T-cell by concanavalin-A. As a cell surface receptor, TREML2 may play a role in the innate and adaptive immune response. TREML2 also acts as a counter-receptor for CD276 and interaction with CD276 on T-cells enhances T-cell activation. It has shown that TREML2 may be involved in the innate immune response based on its expression profile and the fact that it is up-regulated in response to inflammation.

配方:PBS

分子量:54.2 kDa

序列:Gly19-Ser268

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.