



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Leucine-rich alpha-2-glycoprotein Protein

产品货号: TP05066

第三版

别名:Leucine-rich alpha-2-glycoprotein;HMFT1766; LRG; LRG1

描述:Leucine-rich alpha-2-glycoprotein is a secreted protein and contains 8 LRR (leucine-rich) repeats and 1 LRRCT domain. The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation. Levels of the LRG protein are markedly elevated in acute appendicitis and therefore could be used as a diagnostic aid.

配方:PBS

分子量:35.3 kDa

序列:Val36-Gln347

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.