

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Neuritin Protein

产品货号: TP05089

第三版

别名:Neuritin;NRN1;NRN

描述:Neuritin/NRN1 is a member of the neuritin family and can be expressed in postmitotic-differentiating neurons of the developmental nervous system and neuronal structures associated with plasticity in the adult. Neuritin/NRN1 promotes neurite outgrowth, arborization and neuritogenesis. The protein contains a consensus cleavage signal found in glycosylphosphatidylinositol (GPI)-anchored proteins. Overexpression of the encoded protein may be associated with astrocytoma progression.

配方:PBS

分子量:12.1 kDa

序列:Ala28-Gly116

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.