

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PI3K-Interacting Protein 1/PIK3IP1 Protein(C-6His)

产品货号: TP05090

第三版

别名:Kringle domain-containing protein HGFL; PIK3IP1; HGFL

描述:Phosphoinositide-3-kinase-interacting protein 1(PIK3IP1) is an enzyme that in humans is encoded by the PIK3IP1 gene.It is a negative regulator of phosphatidylinositol-3-kinase (PI3K), suppresses the development of hepatocellular carcinoma. The gene encoding PIK3IP1 maps to human chromosome 22, which houses over 500 genes and is the second smallest human chromosome. Mutations in several of the genes that map to chromosome 22 are involved in the development of Phelan-McDermid syndrome, Neurofibromatosis type 2, autism and schizophrenia.

配方:PBS

分子量:16.7 kDa

序列:Ser22-Thr168

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.